


● Puncture Closure Device

- Percutaneous Suturing
- Ligating Abdominal Wall Bleeders
- Trocar Closure without Extending the Incision
- Tacking up Hernia Mesh Prior to Final Placement



- 14-gauge outside diameter.
- Enables the surgeon to rapidly close each trocar site securely and safely.
- Simple and efficient one-hand operation.
- Accommodates any size suture, any size trocar.
- Cost-effective design permits multiple trocar site closure with a single device per case.
- Disposable, 10 units per box.

Catalog #RSG-14

Part #76001-01M

Ordering Instructions



To place an order in Switzerland or to receive more information contact:
Balmer Medical SA
T. +41 24 436 29 07
info@balmermedical.ch



ISO 9001

Precision instruments for specialized applications

Ranfac Corp. - 30 Doherty Ave. - Avon, MA 02322 - USA
US: 1-800-272-6322 - Int'l: 1-508-588-4400

Distributed in Switzerland by:

 **balmer medical**

Balmer Medical SA - Route de Provence 52 - CH-1426 Concise - Switzerland

Tel +41 24 436 29 07 • Fax +41 24 436 29 08 • www.balmermedical.ch • Email: info@balmermedical.ch